




**ZOLL**



**AED**

# ZOLL AED 3™ BLS

## Real CPR Help® with Programme Management Onboard™ Even Better Support for Rescuers

CPR Dashboard™ uses Real CPR Help to show rate and depth of chest compressions in adult as well as child rescues (shown here)—an industry first for an AED.

Easily reconfigure via interactive touchscreen to match clinical protocols and turn audio recording on/off.

Shock button arms within 8 seconds of stopping CPR (with fresh battery).



Programme Management Onboard is cloud-connected to PlusTrac™ for automatic asset tracking log, email notification of self-test failures, as well as pads and battery expiry countdown.

CPR countdown timer shows time remaining in current CPR cycle.

Easily switch to Child mode with the press of a button. No pads to change with CPR Uni-padz™.

To learn more about the ZOLL AED 3™, call us at 01928 595 160  
or visit us at [www.zoll.com/uk/aed3](http://www.zoll.com/uk/aed3)

# ZOLL®

AED 3 is not available for sale in the U.S. or Canada. The product has not received regulatory clearance/approval from the Food and Drug Administration or Health Canada.

©2016 ZOLL Medical Corporation. All rights reserved. CPR Dashboard, PlusTrac, Programme Management Onboard, Real CPR Help, Uni-padz, ZOLL AED 3, and ZOLL are trademarks or registered trademarks of ZOLL Medical Corporation in the United States and/or other countries.

Printed in USA.  
MCN PP 1409 0073-05