

AED Plus® Pelican Case

Technical Specifications



Total Protection for the Best Support




Providing an excellent shock-resistant carrier for your ZOLL® AED Plus®, the Pelican Case and its molded high-density foam provide total protection from the rigors of impact, vibration, and shock. Made of high-impact structural copolymer, it is extremely strong and durable. With its exclusive 0.25" (6.4 mm) neoprene O-ring, and ABS latches, this case seals perfectly. The Pelican Case also includes an automatic purge valve for quick equalization after changes in atmospheric pressure. The AED Plus Pelican Case is NATO-codified and tested to MIL C-4150J (Military Standard) and IP67 (water and dust ingress protection). You break it, we replace the case ... forever.



AED Plus Pelican Case

- Waterproof
- Crushproof
- Dustproof
- Corrosion proof
- Chemical resistant

ZOLL®



ZOLL's AED Plus Large Pelican Case allows you to store your AED Plus and an extra set of CPR-D-padz with the option of also storing a set of Pedi-padz II.

**ADVANCING
RESUSCITATION.
TODAY.®**

ZOLL Medical Corporation
Worldwide Headquarters
269 Mill Road
Chelmsford, MA 01824
978-421-9655
800-804-4356

For subsidiary addresses and fax numbers, as well as other global locations, please go to www.zoll.com/contacts.

AED Plus Pelican Case Specifications

ZOLL offers two cases to meet your needs: a large case for easy storage of CPR-D-padz® and Pedi-padz®II or a smaller case for the AED Plus alone. Each case comes complete with custom-molded foam inserts for a perfect fit. When storing an AED Plus in a pelican case, the status indicator is not visible and the audible alert sound may be muffled. ZOLL recommends opening the case at least once a month to confirm device status.

Large Case (8000-0837-01)

Large Case with cutouts for AED Plus, CPR-D-padz, or Pedi-padz II

Inside Dimensions: 17.9" x 12.8" x 6.8"
(45.5 cm x 31.8 cm x 17.3 cm)

Outside Dimensions: 19" x 15.4" x 7.6"
(48.3 cm x 39.1 cm x 19.3 cm)

Temperature Ratings: Min: -10° F
(-23° C) Max: 210° F (99° C)

Buoyancy: Floats in salt water with 40 lbs. (18.1 kg) load

Color: Black

Approvals: IP67 Military Standard
C-4150J

Small (8000-0836-01)

Small Case with cutouts for AED Plus only

Inside Dimensions: 14.7" x 10.4" x 6"
(37.3 cm x 26.4 cm x 15.2 cm)

Outside Dimensions: 16" x 13" x 6.9"
(40.6 cm x 33 cm x 17.5 cm)

Temperature Ratings: Min: -10° F (-23° C)
Max: 210° F (99° C)

Buoyancy: Floats in salt water with 20 lbs. (9.07 kg) load

Color: Black

Approvals: IP67 Military Standard
C-4150J

A Rugged Case for the First and Only Full-Rescue AED



Specifications subject to change without notice.

©2015 ZOLL Medical Corporation. All rights reserved. Advancing Resuscitation Today, AED Plus, CPR-D-padz, Pedi-padz II, and ZOLL are trademarks or registered trademarks of ZOLL Medical Corporation in the United States and/or other countries. All other trademarks are the property of their respective owners.

Printed in U.S.A.

MCN PP 1503 0115

ZOLL®